

CHAPRO

Compression Hearing-Aid Processing Library API Documentation

**Boys Town National Research Hospital
11 June 2018**

TABLE OF CONTENTS

CHAPRO LIBRARY OVERVIEW	3
CHAPRO FUNCTION DESCRIPTIONS	6
cha_allocate.....	6
cha_cleanup.....	7
cha_data_gen.....	8
cha_hex_patch.....	9
cha_scale	10
cha_fft_cr	11
cha_fft_rc	12
cha_fft	13
cha_ifft	14
cha_prepare	15
cha_version	16
cha_agc_prepare	17
cha_agc_input	18
cha_agc_channel	19
cha_agc_output	20
cha_icmp_prepare	21
cha_icmp_process	22
cha_afc_prepare	23
cha_afc_input.....	24
cha_afc_output.....	25
cha_ciirfb_design.....	26
cha_ciirfb_prepare	27
cha_ciirfb_analyze	28
cha_ciirfb_synthesize.....	29
cha_cfirfb_prepare	30
cha_cfirfb_analyze.....	31
cha_cfirfb_synthesize	32
cha_firfb_prepare.....	33
cha_firfb_analyze.....	34
cha_firfb_synthesize	35
cha_iirfb_design.....	36
cha_iirfb_prepare	37
cha_iirfb_analyze.....	38
cha_firfb_synthesize	39
Appendix A. Test programs	40
Appendix B. CLS Prescription	41
Appendix C. DSL Prescription	42
Appendix D. WDRC Parameters	43

CHAPRO LIBRARY OVERVIEW

CHAPRO is a library of functions that may be used to implement simulations of compression hearing-aid signal processing. Four different types of signal processing strategies are included:

- (1) FIR filter-bank frequency analysis with automatic gain control;
- (2) IIR filter-bank frequency analysis with automatic gain control;
- (3) Complex FIR filter-bank frequency analysis with instantaneous compression;
- (4) Complex IIR filter-bank frequency analysis with instantaneous compression.

A modular design has been adopted to facilitate replacement of library functions with alternative signal-processing implementations. Each of the six major modules contains (1) a *preparation* function that allocates memory and initializes variables and (2) one or more *processing* functions that perform signal processing.

1. FIR filter-bank
 - a. `firfb_prepare`
 - b. `firfb_analyze`
 - c. `firfb_synthesize`
2. IIR filter-bank
 - a. `iirfb_design`
 - b. `iirfb_prepare`
 - c. `iirfb_analyze`
 - d. `iirfb_synthesize`
3. Automatic gain control
 - a. `agc_prepare`
 - b. `agc_input`
 - c. `agc_channel`
 - d. `agc_output`
4. Adaptive feedback cancelation
 - a. `afc_prepare`
 - b. `afc_input`
 - c. `afc_output`
5. Complex IIR filter-bank
 - a. `ciirfb_design`
 - b. `ciirfb_prepare`
 - c. `ciirfb_analyze`
 - d. `ciirfb_synthesize`
6. Instantaneous compression
 - a. `icmp_prepare`
 - b. `icmp_process`

For each module, variables are initialized and associated data memory is allocated by the preparation functions. This storage is combined into a single data structure to facilitate creation of firmware for real-time implementation on signal-processing hardware that might not have an operating system. The CHAPRO library includes a function that generates a C-code representation of this initialized data.

- `data_gen`

For desktop simulation, the CHAPRO library includes core functions for memory allocation and disposal.

- `prepare`
- `allocate`
- `cleanup`

Functions for FFT of real signals are included among the core functions.

- `fft_rc`
- `fft_cr`

Finally, the CHAPRO library includes a function that returns a version description string.

- `version`

To simulate complex IIR filter-bank frequency analysis with instantaneous compression, variable initialization and memory allocation is performed by calling the following functions.

- `ciirfb_prepare`
- `icmp_prepare`

Subsequent signal processing is performed by calling the following functions.

- `ciirfb_analyze`
- `icmp_process`
- `ciirfb_synthesize`

Several examples test basic aspects of these functions.

- `tst_cifa` – tests complex IIR filter-bank analysis
- `tst_cifio` – tests simple waveform complex IIR processing
- `tst_cifsc` – tests simple waveform complex IIR processing with soundcard

To simulate FIR filter-bank frequency analysis with AGC compression, variable initialization and memory allocation is performed by calling the following prepare functions.

- `firfb_prepare`
- `agc_prepare`

Subsequent signal processing is performed by calling the following process functions.

- `agc_input`
- `firfb_analyze`
- `agc_channel`
- `firfb_synthesize`

- `agc_output`

Several examples test basic aspects of these functions.

- `tst_ffa` – tests filter-bank analysis
- `tst_ffio` – tests simple waveform FIR processing
- `tst_ffsc` – tests speech-waveform FIR & AGC with soundcard

To simulate IIR filter-bank frequency analysis with AGC compression and adaptive feedback cancelation, variable initialization and memory allocation is performed by calling the following functions.

- `afc_prepare`
- `agc_prepare`
- `iirfb_prepare`

Subsequent signal processing is performed by calling the following functions.

- `afc_input`
- `agc_input`
- `iirfb_analyze`
- `agc_channel`
- `iirfb_synthesize`
- `agc_output`
- `afc_output`

Several examples test basic aspects of these functions.

- `tst_ifa` – tests filter-bank analysis
- `tst_ifio` – tests simple waveform IIR processing
- `tst_ifsc` – tests speech-waveform IIR & AGC with soundcard
- `tst_iffb` – tests speech-waveform IIR & AFC with soundcard
- `tst_gha` – tests speech-waveform IIR & AFC & AGC with soundcard

All examples require the SIGPRO library from BTNRH (<http://audres.org/rc/sigpro>). The soundcard examples also require the ARSC library (<http://audres.org/rc/arsc>).

CHAPRO FUNCTION DESCRIPTIONS

cha_allocate

Allocates memory attached to CHAPRO data structure.

(void) **cha_allocate**(CHA_PTR **cp**, int **cnt**, int **siz**, int **idx**)

Function arguments

cp	pointer to CHAPRO data structure
cnt	Number of elements to allocate.
siz	Size of each element.
idx	Index into CHAPPRO data structure.

Return Value

none

Remarks

A pointer to the allocated memory is stored in the CHAPRO data structure at the location specified by **idx**.

See Also

cha_cleanup

cha_cleanup

Frees all memory attached to CHAPRO data structure.

(void) **cha_cleanup**(CHA_PTR **cp**)

Function arguments

cp pointer to CHAPRO data structure

Return Value

none

Remarks

Should always be the last function called in the CHAPRO library.

See Also

cha_allocate

cha_data_gen

Generates C code that represents the CHAPRO data structure.

(int) **cha_data_gen**(CHA_PTR **cp**, char ***fn**)

Function arguments

cp	pointer to CHAPRO data structure
fn	Pointer to output filename.

Return Value

Error code:

- 0 – no error
- 1 – can't open output file
- 2 – data structure not yet initialized
- 3 – data structure contains no data

Remarks

The C code generated by this function represents the CHAPRO data structure after variables have been initialized and data memory has been allocated by prior calls to any preparation functions. The code is written to the file specified by **fn**.

cha_hex_patch

Patch hex file by replacing cha_data.

(int) **cha_hex_patch**(CHA_PTR **cp**, char ***ifn**, char ***ofn**)

Function arguments

cp	pointer to CHAPRO data structure
ifn	Pointer to input filename.
ofn	Pointer to output filename.

Return Value

Error code:
0 – no error
1 – can't open output file
2 – data structure not yet initialized
3 – data structure contains no data

Remarks

Copies the input hex file **ofn** to an output hex file **ifn** with replacement of cha_data.

cha_scale

Applies scales factor to a chunk of the input or output stream.

(void) **cha_scale**(float ***x**, int **cs**, float **scale**)

Function arguments

x	pointer to input signal
cs	chunk size
scale	pointer to CHAPRO data structure

Return Value

None.

Remarks

The scaled output signal overwrites the input signal.

cha_fft_cr

Inverse Fourier transform complex frequency components into real signal.

(void) **cha_fft_cr**(float ***x**, int **n**)

Function arguments

x	Complex frequency components are replaced by real-valued signal.
n	Number of points in the signal.

Return Value

None

Remarks

The input array must be dimensioned to accommodate **n+2** float values. The number of complex frequency components is **(n+2)/2**.

cha_fft_rc

Fourier transform real signal into complex frequency components.

(void) **cha_fft_rc**(float ***x**, int **n**)

Function arguments

x	Real-valued signal is replaced by complex frequency components
n	Number of points in the signal.

Return Value

None

Remarks

The input array must be dimensioned to accommodate **n+2** float values. The number of complex frequency components is **(n+2)/2**.

cha_fft

Fourier transform a complex time signal to complex frequency components.

(void) **cha_fft**(float ***x**, int **n**)

Function arguments

x	Complex frequency components are replaced by a complex signal.
n	Number of points in the signal.

Return Value

None

Remarks

The input array must be dimensioned to accommodate **n**×**2** float values.

cha_ift

Inverse Fourier transform complex frequency components into a complex-valued signal.

(void) **cha_ift**(float ***x**, int **n**)

Function arguments

x	Complex signal is replaced by complex frequency components
n	Number of points in the signal.

Return Value

None

Remarks

The input array must be dimensioned to accommodate **n**×2 float values.

cha_prepare

CHAPRO data structure preparation function.

(void) **cha_prepare**(CHA_PTR **cp**)

Function arguments

cp pointer to CHAPRO data structure

Return Value

None.

Remarks

Should be called only once and prior to calling other library functions; however, violations of this rule may be tolerated.

cha_version

Returns a string that describes the current version of the CHAPRO library.

(char *) **cha_version**(void)

Function arguments

none

Return Value

Pointer to version string.

Remarks

An example of the return value, "CHAPro version 0.03, 6-Nov-2016".

cha_agc_prepare

Automatic-gain-control preparation function.

(int) **cha_agc_prepare**(CHA_PTR **cp**, CHA_DSL ***dsl**, CHA_WDRC ***gha**)

Function arguments

cp	pointer to CHA data structure
dsl	pointer to DSL prescription structure (see Appendix C)
gha	pointer to WDRC prescription structure (see Appendix D)

Return Value

Error code:
0 – no error

Remarks

Initializes variables and allocates memory for automatic gain control. Chunk size is the number of samples read from the input signal and written to the output signal with each call to `cha_agc_process`.

See Also

cha_agc_process

cha_agc_input

Automatic-gain-control processing function.

(void) **cha_agc_input**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs single-channel, automatic-gain-control processing on CHAPRO input signal. Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_agc_prepare, cha_agc_output

cha_agc_channel

Automatic-gain-control processing function.

(void) **cha_agc_channel**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs multi-channel, automatic-gain-control processing. Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_agc_prepare

cha_agc_output

Automatic-gain-control processing function.

(void) **cha_agc_output**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs single-channel, automatic-gain-control processing on CHAPRO output signal. Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_agc_prepare, **cha_agc_input**

cha_icmp_prepare

Instantaneous-compression preparation function.

(int) **cha_icmp_prepare**(CHA_PTR **cp**,
float ***Lc**, float ***Gc**, double **lr**, int **np**, int **ds**)

Function arguments

cp	pointer to CHA data structure
Lc	pointer to level array
Lc	pointer to gain array
lr	level reference
np	number of points in level and gain arrays
ds	down-sample factor

Return Value

Error code:
0 – no error

Remarks

Initializes variables and allocates memory for instantaneous compression. Chunk size is the number of samples read from the input signal and written to the output signal with each call to `cha_icmp_process`.

See Also

cha_icmp_process

cha_icmp_process

Instantaneous-compression processing function.

(void) **cha_icmp_process**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs instantaneous compression. Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_icmp_prepare

cha_afc_prepare

Configure feedback management.

(int) **cha_afc_prepare**(CHA_PTR **cp**, double **mu**, double **rho**, int **afl**,
int **wfl**, int **pfl**, int **hdel**, double **fbg**, int **sqm**)

Function arguments

cp	pointer to CHAPRO data structure
mu	AFC filter-estimation step size
rho	AFC filter-estimation forgetting factor
eps	AFC filter-estimation power threshold
afl	length of adaptive-feedback-filter response
wfl	length of signal-whitening-filter response
pfl	length of persistent-feedback-filter response
hdel	output-to-input hardware delay (samples)
fbg	simulated-feedback gain
sqm	option to save quality metric

Return Value

Error code:
0 – no error

Remarks

Initializes variables and allocates memory for adaptive feedback cancelation (AFC), which is controlled by parameters **mu**, **rho**, **eps**, and **afl**. An optional signal-whitening filter is specified by setting a positive value for **wfl**. An optional persistent-feedback filter is specified by setting a positive value for **pfl**. (These optional filters are not yet implemented.) For testing purposes, simulated feedback is enabled by setting **fbg**=1 or disabled by setting **fbg**=0. When simulation is enabled, the feedback-filter misalignment error is saved as a quality metric by setting **sqm**=1.

See Also

cha_afc_input, **cha_afc_output**

cha_afc_input

Process input signal to remove feedback.

(void) **cha_afc_input**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHAPRO data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Removes estimated feedback from input signal. Chunk size is the number of samples read from the input signal and written to the output signal. Optionally simulates feedback and saves misalignment error as a quality metric.

See Also

cha_afc_prepare, cha_afc_output

cha_afc_output

Save output signal for feedback management.

(void) **cha_afc_output**(CHA_PTR **cp**, float ***x**, int **cs**)

Function arguments

cp	pointer to CHAPRO data structure
x	pointer to input signal
cs	chunk size

Return Value

none

Remarks

Assists feedback management by saving the output of the hearing-aid processing. Chunk size is the number of samples read from the input signal.

See Also

cha_afc_prepare, cha_afc_input

cha_ciirfb_design

Complex IIR filter-bank design function.

```
(int) cha_ciirfb_design(float *z, float *p, float *g, int *d, int nc,
                        double *fc, double *bw, double sr, double td)
```

Function arguments

z	pointer to IIR filter complex zeros
p	pointer to IIR filter complex poles
g	pointer IIR filter complex gain
d	pointer IIR filter delay
nc	number of frequency bands
fc	pointer to list of center frequencies (Hz)
bw	pointer to list of bandwidths (Hz)
sr	sampling rate (samples/second)
td	target group delay (ms)

Return Value

Error code:
0 – no error

Remarks

Computes zeros and poles for the complex IIR filter-bank. The filterbank design is based on fourth-order gammatone bandpass filters. The center frequency (**fc**) and bandwidth (**bw**) arrays are input arguments and their size is equal to the number of frequency bands (**nc**). The zero (**z**) and pole (**p**) arrays are output arguments and are complex numbers, so are stored as sequential real and imaginary parts. The size of the zeros & poles arrays is two times the product of the number of frequency bands (**nc**) and the number of zeros (and poles) per band, which is always equal to 4. The gain array (**g**) is complex, so its size two times the number of frequency bands (**nc**). The size of the filter delay array (**d**) is equal to the number of frequency bands (**nc**). The filter-design process attempts to align the impulse response of each frequency band to the target delay (**td**).

See Also

cha_ciirfb_prepare, **cha_ciirfb_analyze**, **cha_ciirfb_synthesize**

cha_ciirfb_prepare

Complex IIR filter-bank preparation function.

```
(int) cha_ciirfb_prepare(CHA_PTR cp, float *z, float *p, float *g, int *d,
    int nc, int nz, double sr, int cs)
```

Function arguments

cp	pointer to CHA data structure
z	pointer to IIR filter complex zeros
p	pointer to IIR filter complex poles
g	pointer IIR filter complex gain
d	pointer IIR filter delay
nc	number of frequency bands
nz	number of zeros (and poles) for each band
sr	sampling rate (samples/second)
cs	chunk size

Return Value

Error code:
0 – no error

Remarks

Initializes variables and allocates memory for the IIR filter-bank. The zero (**z**) and pole (**p**) arrays are output arguments and are complex numbers, so are stored as sequential real and imaginary parts. The size of the zero & pole arrays is two times the product of the number of frequency bands (**nc**) and the number of zeros per band (**nz**). The number of zeros (and poles) per band (**nz**) should be an even number and should include conjugate pairs, so that second-order sections will have real coefficients. The gain array (**g**) is complex, so its size two times the number of frequency bands (**nc**). The size of the filter delay array (**d**) is equal to the number of frequency bands (**nc**). Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_ciirfb_design, **cha_ciirfb_analyze**, **cha_ciirfb_synthesize**

cha_ciirfb_analyze

Complex IIR filter-bank frequency-analysis function.

(void) **cha_ciirfb_analyze**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to complex output signal
cs	chunk size

Return Value

none

Remarks

Performs complex IIR filter-bank analysis. Chunk size is the number of samples read from the input and written to the output.

See Also

cha_ciirfb_prepare, cha_ciirfb_synthesize

cha_ciirfb_synthesize

Complex IIR filter-bank frequency-synthesis function.

(void) **cha_ciirfb_synthesize**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to complex input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs complex filter-bank frequency-synthesis. Frequency bands are summed and only the real part is output. Chunk size is the number of samples read from the input and written to the output.

See Also

cha_ciirfb_prepare, **cha_ciirfb_analyze**

cha_cfirfb_prepare

Complex FIR filter-bank preparation function.

(int) **cha_cfirfb_prepare**(CHA_PTR **cp**, double ***cf**, int **nc**, double **sr**,
int **nw**, int **wt**, int **cs**)

Function arguments

cp	pointer to CHA data structure
cf	list frequency band edges (kHz)
nc	number of frequency bands
sr	sampling rate (samples/second)
nw	window size (samples)
wt	window type (0=Hamming, 1=Blackman)
cs	chunk size

Return Value

Error code:
0 – no error

Remarks

Initializes variables and allocates memory for the complex FIR filter-bank. Chunk size is the number of samples read from the input and written to the output.

See Also

cha_cfirfb_analyze, **cha_cfirfb_synthesize**

cha_cfirfb_analyze

Complex FIR filter-bank frequency-analysis function.

(void) **cha_cfirfb_analyze**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to real input signal
y	pointer to complex output signal
cs	chunk size

Return Value

none

Remarks

Performs complex FIR filter-bank analysis. Chunk size is the number of samples read from the input and written to the output.

See Also

cha_cfirfb_prepare, **cha_cfirfb_synthesize**

cha_cfirfb_synthesize

Complex FIR filter-bank frequency-synthesis function.

(void) **cha_cfirfb_synthesize**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to complex input signal
y	pointer to real output signal
cs	chunk size

Return Value

none

Remarks

Performs complex FIR filter-bank synthesis. Chunk size is the number of samples read from the input and written to the output.

See Also

cha_cfirfb_prepare, **cha_cfirfb_analyze**

cha_firfb_prepare

FIR filter-bank preparation function.

(int) **cha_firfb_prepare**(CHA_PTR **cp**,
double ***cf**, int **nc**, double **sr**, int **nw**, int **wt**, int **cs**)

Function arguments

cp	pointer to CHA data structure
cf	list frequency band edges (kHz)
nc	number of frequency bands
sr	sampling rate (samples/second)
nw	window size (samples)
wt	window type (0=Hamming, 1=Blackman)
cs	chunk size

Return Value

Error code:
0 – no error

Remarks

Initializes variables and allocates memory for the FIR filter-bank. Chunk size is the number of samples read from the input and written to the output.

See Also

cha_firfb_analyze, cha_firfb_synthesize

cha_firfb_analyze

FIR filter-bank frequency-analysis function.

(void) **cha_firfb_analyze**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs FIR filter-bank analysis. Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_firfb_prepare, **cha_firfb_synthesize**

cha_firfb_synthesize

FIR filter-bank frequency-synthesis function.

(void) **cha_firfb_synthesize**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs FIR filter-bank synthesis. Chunk size is the number of samples read from the input signal and written to the output signal.

cha_iirfb_design

IIR filter-bank design function.

```
(int) cha_iirfb_design(float *z, float *p, float *g, int *d,
    int nc, int nz, double sr, double td)
```

Function arguments

z	pointer to IIR filter complex zeros
p	pointer to IIR filter complex poles
g	pointer to IIR filter gain
d	pointer to IIR filter delay
cf	pointer to cross-over frequencies
nc	number of frequency bands
nz	number of zeros (and poles) for each band
sr	sampling rate (samples/second)
td	impulse response target delay (millisecond)

Return Value

Error code:
0 – no error

Remarks

Computes zeros and poles for the IIR filter-bank. The filterbank design is based on Butterworth bandpass filters sandwiched between Butterworth low-pass and high-pass filters. The array of cross-over frequencies (**cf**) is an input argument and its size is one less than the number of frequency bands (**nc**). The zero (**z**) and pole (**p**) arrays are output arguments and are complex numbers, so are stored as sequential real and imaginary parts. The size of the zeros & poles arrays is two times the product of the number of frequency bands (**nc**) and the number of zeros per band (**nz**). The number of zeros (and poles) per band (**nz**) should be an even number and should include conjugate pairs, so that second-order sections have real coefficients. The number of filter gains and filter delays (**g** & **d**) is the number of frequency bands (**nc**). The filter-design process attempts to align the impulse response of each frequency band to the target delay (**td**).

See Also

cha_iirfb_prepare, **cha_iirfb_analyze**, **cha_iirfb_synthesize**

cha_iirfb_prepare

IIR filter-bank preparation function.

```
(int) cha_iirfb_prepare(CHA_PTR cp, float *z, float *p, float *g, int *d,
    int nc, int nz, double sr, int cs)
```

Function arguments

cp	pointer to CHA data structure
z	pointer to IIR filter complex zeros
p	pointer to IIR filter complex poles
g	pointer IIR filter gain
d	pointer IIR filter delay
nc	number of frequency bands
nz	number of zeros (and poles) for each band
sr	sampling rate (samples/second)
cs	chunk size

Return Value

Error code:
0 – no error

Remarks

Initializes variables and allocates memory for the IIR filter-bank. The zero (**z**) and pole (**p**) arrays are output arguments and are complex numbers, so are stored as sequential real and imaginary parts. The size of the zero & pole arrays is two times the product of the number of frequency bands (**nc**) and the number of zeros per band (**nz**). The number of zeros (and poles) per band (**nz**) should be an even number and should include conjugate pairs, so that second-order sections will have real coefficients. The filter gains and filter delays (**g** & **d**) are input arguments and have array sizes equal to the number of frequency bands (**nc**). Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_iirfb_design, **cha_iirfb_analyze**, **cha_iirfb_synthesize**

cha_iirfb_analyze

IIR filter-bank frequency-analysis function.

(void) **cha_iirfb_analyze**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs IIR filter-bank analysis. Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_iirfb_prepare, **cha_iirfb_synthesize**

cha_firfb_synthesize

IIR filter-bank frequency-synthesis function.

(void) **cha_iirfb_synthesize**(CHA_PTR **cp**, float ***x**, float ***y**, int **cs**)

Function arguments

cp	pointer to CHA data structure
x	pointer to input signal
y	pointer to output signal
cs	chunk size

Return Value

none

Remarks

Performs IIR filter-bank synthesis. Chunk size is the number of samples read from the input signal and written to the output signal.

See Also

cha_iirfb_prepare, cha_iirfb_analyze

Appendix A. Test programs

Several examples that test basic aspects of complex IIR filter-bank and instantaneous-compression.

- `tst_cifa` – tests complex IIR filter-bank analysis
- `tst_cifio` – tests simple waveform complex IIR processing
- `tst_cifsc` – tests simple waveform complex IIR processing with soundcard

Several examples that test basic aspects of FIR filter-bank and automatic-gain-control.

- `tst_ffa` – tests FIR filter-bank analysis
- `tst_ffio` – tests simple waveform FIR processing
- `tst_ffsc` – tests speech-waveform FIR & AGC with soundcard

Several examples that test basic aspects of IIR filter-bank with automatic-gain-control and adaptive feedback cancelation.

- `tst_ifa` – tests IIR filter-bank analysis
- `tst_ifio` – tests simple waveform IIR processing
- `tst_ifsc` – tests speech-waveform IIR & AGC
- `tst_iffb` – tests speech-waveform IIR & AFC with soundcard
- `tst_gha` – tests speech-waveform IIR & AFC & AGC with soundcard

These test programs are all written in C and produce results that are written to a subfolder called “test.” Each test program has a corresponding MATLAB script for viewing the results. The test programs that contain IIR processing require filter coefficients to be precomputed by a MATLAB script called `iirfb`.

Appendix B. CLS Prescription

Structure CHA_CLS specifies the CLS prescription.

```
#define CLS_MXCH 32          // maximum number of channels

typedef struct {
    int cm;                  // compression mode
    int nc;                  // number of channels
    double fc[CLS_MXCH];    // center frequency
    double bw[CLS_MXCH];    // bandwidth
    double Gcs[CLS_MXCH];   // gain at compression start
    double Gcm[CLS_MXCH];   // gain at compression middle
    double Gce[CLS_MXCH];   // gain at compression end
    double Gmx[CLS_MXCH];   // maximum gain
    double Lcs[CLS_MXCH];   // level at compression start
    double Lcm[CLS_MXCH];   // level at compression middle
    double Lce[CLS_MXCH];   // level at compression end
    double Lmx[CLS_MXCH];   // maximum output level
} CHA_CLS;
```

Appendix C. DSL Prescription

Structure `CHA_DSL` specifies the DSL prescription.

```
#define DSL_MXCH 32                // maximum number of channels

typedef struct {
    double attack;                  // attack time (ms)
    double release;                 // release time (ms)
    double maxdB;                   // maximum output (dB SPL)
    int ear;                        // 0=left, 1=right
    int nchannel;                   // number of channels
    double cross_freq[DSL_MXCH];   // cross frequencies (Hz)
    double tkgain[DSL_MXCH];       // compression-start gain
    double cr[DSL_MXCH];           // compression ratio
    double tk[DSL_MXCH];           // compression-start kneepoint
    double bolt[DSL_MXCH];         // broadband output limiting threshold
} CHA_DSL;
```

Appendix D. WDRC Parameters

Structure `CHA_WDRC` specifies single-channel WDRC parameters

```
typedef struct {  
    double attack;           // attack time (ms)  
    double release;          // release time (ms)  
    double fs;               // sampling rate (Hz)  
    double maxdB;            // maximum signal (dB SPL)  
    double tkgain;           // compression-start gain  
    double tk;               // compression-start kneepoint  
    double cr;               // compression ratio  
    double bolt;             // broadband output limiting threshold  
} CHA_WDRC;
```