

# **CHAPRO**

## **Compression Hearing-Aid Processing Library API Documentation**

**Boys Town National Research Hospital  
27 April 2018**

**TABLE OF CONTENTS**

CHAPRO LIBRARY OVERVIEW .....	3
CHAPRO FUNCTION DESCRIPTIONS .....	6
cha_allocate.....	6
cha_cleanup.....	7
cha_data_gen.....	8
cha_hex_patch.....	9
cha_scale .....	10
cha_fft_cr .....	11
cha_fft_rc .....	12
cha_fft .....	13
cha_ifft .....	14
cha_prepare .....	15
cha_version .....	16
cha_agc_prepare .....	17
cha_agc_input .....	18
cha_agc_channel .....	19
cha_agc_output .....	20
cha_compressor_prepare.....	21
cha_compressor_process.....	22
cha_afc_prepare .....	23
cha_afc_input.....	24
cha_afc_output.....	25
cha_cgtfb_prepare.....	26
cha_cgtfb_analyze.....	27
cha_cgtfb_synthesize .....	28
cha_cfirfb_prepare .....	29
cha_cfirfb_analyze.....	30
cha_cfirfb_synthesize .....	31
cha_firfb_prepare .....	32
cha_firfb_analyze.....	33
cha_firfb_synthesize .....	34
cha_iirfb_prepare .....	35
cha_iirfb_analyze.....	36
cha_firfb_synthesize .....	37
Appendix A. Test programs.....	38
Appendix B. CLS Prescription .....	39
Appendix C. DSL Prescription .....	40
Appendix D. WDRC Parameters .....	41

## CHAPRO LIBRARY OVERVIEW

CHAPRO is a library of functions that may be used to implement simulations of compression hearing-aid signal processing. Two different types of signal processing strategies are included: (1) gammatone filter-bank frequency analysis with instantaneous compression and (2) FIR filter-bank frequency analysis with automatic gain control.

A modular design has been adopted to facilitate replacement of library functions with alternative signal-processing implementations. Each of the five major modules contains (1) a *preparation* function that allocates memory and initializes variables and (2) one or more *processing* functions that perform signal processing.

1. Complex Gammatone filter-bank
  - a. `cgtfb_prepare`
  - b. `cgtfb_analyze`
  - c. `cgtfb_synthesize`
2. Instantaneous compression
  - a. `compressor_prepare`
  - b. `compressor_process`
3. FIR filter-bank Fourier Transform
  - a. `firfb_prepare`
  - b. `firfb_analyze`
  - c. `firfb_synthesize`
4. IIR filter-bank Fourier Transform
  - a. `iirfb_prepare`
  - b. `iirfb_analyze`
  - c. `iirfb_synthesize`
5. Automatic gain control
  - a. `agc_prepare`
  - b. `agc_input`
  - c. `agc_channel`
  - d. `agc_output`
6. Adaptive feedback cancelation
  - a. `afc_prepare`
  - b. `afc_input`
  - c. `afc_output`

All variables are initialized and related data memory is allocated by the preparation functions. This storage is combined into a single data structure to facilitate the creation of firmware for real-time implementation on signal-processing hardware. The CHAPRO library includes a function that generates a C-code representation of this initialized data.

- `data_gen`

For desktop simulation, the CHAPRO library includes core functions for memory allocation and disposal.

- `prepare`

- `allocate`
- `cleanup`

Functions for FFT of real signals are included among the core functions.

- `fft_rc`
- `fft_cr`

Finally, the CHAPRO library includes a function that returns a version description string.

- `version`

To simulate gammatone filter-bank frequency analysis with instantaneous compression, variable initialization and memory allocation is performed by calling the following functions.

- `cgtfb_prepare`
- `compressor_prepare`

Subsequent signal processing is performed by calling the following functions.

- `cgtfb_analyze`
- `compressor_process`
- `cgtfb_synthesize`

Several examples test basic aspects of these functions.

- `tst_gfa` – tests gammatone filter-bank analysis
- `tst_gfio` – tests simple waveform gammatone processing
- `tst_gfsc` – tests simple waveform gammatone processing with soundcard

To simulate FIR filter-bank frequency analysis with AGC compression, variable initialization and memory allocation is performed by calling the following functions.

- `firfb_prepare`
- `agc_prepare`

Subsequent signal processing is performed by calling the following functions.

- `agc_input`
- `firfb_analyze`
- `agc_channel`
- `firfb_synthesize`
- `agc_output`
- `feedback_record`

Several examples test basic aspects of these functions.

- `tst_ffa` – tests filter-bank analysis
- `tst_ffio` – tests simple waveform FIR processing
- `tst_ffsc` – tests speech-waveform FIR & AGC with soundcard

To simulate IIR filter-bank frequency analysis with AGC compression and adaptive feedback cancelation, variable initialization and memory allocation is performed by calling the following functions.

- `agc_prepare`

- `iirfb_prepare`
- `afc_prepare`

Subsequent signal processing is performed by calling the following functions.

- `afc_input`
- `agc_input`
- `iirfb_analyze`
- `agc_channel`
- `iirfb_synthesize`
- `agc_output`
- `afc_output`

Several examples test basic aspects of these functions.

- `tst_ifa` – tests filter-bank analysis
- `tst_ifio` – tests simple waveform IIR processing
- `tst_ifsc` – tests speech-waveform IIR & AGC with soundcard
- `tst_iffb` – tests speech-waveform IIR & AFC with soundcard
- `tst_gha` – tests speech-waveform IIR & AFC & AGC with soundcard

All examples require the SIGPRO library from BTNRH (<http://audres.org/rc/sigpro>). The soundcard examples also require the ARSC library (<http://audres.org/rc/arsc>).

## CHAPRO FUNCTION DESCRIPTIONS

### *cha\_allocate*

Allocates memory attached to CHAPRO data structure.

(void) **cha\_allocate**(CHA\_PTR **cp**, int **cnt**, int **siz**, int **idx**)

#### Function arguments

<b>cp</b>	pointer to CHAPRO data structure
<b>cnt</b>	Number of elements to allocate.
<b>siz</b>	Size of each element.
<b>idx</b>	Index into CHAPPRO data structure.

#### Return Value

none

#### Remarks

A pointer to the allocated memory is stored in the CHAPRO data structure at the location specified by **idx**.

#### See Also

**cha\_cleanup**

***cha\_cleanup***

Frees all memory attached to CHAPRO data structure.

(void) **cha\_cleanup**(CHA\_PTR **cp**)

**Function arguments**

**cp**                    pointer to CHAPRO data structure

**Return Value**

none

**Remarks**

Should always be the last function called in the CHAPRO library.

**See Also**

**cha\_allocate**

***cha\_data\_gen***

Generates C code that represents the CHAPRO data structure.

(int) **cha\_data\_gen**(CHA\_PTR **cp**, char \***fn**)

**Function arguments**

<b>cp</b>	pointer to CHAPRO data structure
<b>fn</b>	Pointer to output filename.

**Return Value**

Error code:

- 0 – no error
- 1 – can't open output file
- 2 – data structure not yet initialized
- 3 – data structure contains no data

**Remarks**

The C code generated by this function represents the CHAPRO data structure after variables have been initialized and data memory has been allocated by prior calls to any preparation functions. The code is written to the file specified by **fn**.



***cha\_hex\_patch***

Patch hex file by replacing cha\_data.

(int) **cha\_hex\_patch**(CHA\_PTR **cp**, char \***ifn**, char \***ofn**)

**Function arguments**

<b>cp</b>	pointer to CHAPRO data structure
<b>ifn</b>	Pointer to input filename.
<b>ofn</b>	Pointer to output filename.

**Return Value**

Error code:  
0 – no error  
1 – can't open output file  
2 – data structure not yet initialized  
3 – data structure contains no data

**Remarks**

Copies the input hex file **ofn** to an output hex file **ifn** with replacement of cha\_data.

***cha\_scale***

Applies scales factor to a chunk of the input or output stream.

(void) **cha\_scale**(float \***x**, int **cs**, float **scale**)

**Function arguments**

<b>x</b>	pointer to input signal
<b>cs</b>	chunk size
<b>scale</b>	pointer to CHAPRO data structure

**Return Value**

None.

**Remarks**

The scaled output signal overwrites the input signal.

***cha\_fft\_cr***

Inverse Fourier transform complex frequency components into real signal.

(void) **cha\_fft\_cr**(float \***x**, int **n**)

**Function arguments**

<b>x</b>	Complex frequency components are replaced by real-valued signal.
<b>n</b>	Number of points in the signal.

**Return Value**

None

**Remarks**

The input array must be dimensioned to accommodate **n+2** float values. The number of complex frequency components is **(n+2)/2**.

***cha\_fft\_rc***

Fourier transform real signal into complex frequency components.

(void) **cha\_fft\_rc**(float \***x**, int **n**)

**Function arguments**

<b>x</b>	Real-valued signal is replaced by complex frequency components
<b>n</b>	Number of points in the signal.

**Return Value**

None

**Remarks**

The input array must be dimensioned to accommodate **n+2** float values. The number of complex frequency components is **(n+2)/2**.

***cha\_fft***

Fourier transform a complex time signal to complex frequency components.

(void) **cha\_fft**(float \***x**, int **n**)

**Function arguments**

<b>x</b>	Complex frequency components are replaced by a complex signal.
<b>n</b>	Number of points in the signal.

**Return Value**

None

**Remarks**

The input array must be dimensioned to accommodate **n**×**2** float values.

***cha\_ift***

Inverse Fourier transform complex frequency components into a complex-valued signal.

(void) **cha\_ift**(float \***x**, int **n**)

**Function arguments**

<b>x</b>	Complex signal is replaced by complex frequency components
<b>n</b>	Number of points in the signal.

**Return Value**

None

**Remarks**

The input array must be dimensioned to accommodate **n**×2 float values.

***cha\_prepare***

CHAPRO data structure preparation function.

(void) **cha\_prepare**(CHA\_PTR **cp**)

**Function arguments**

**cp**                    pointer to CHAPRO data structure

**Return Value**

None.

**Remarks**

Should be called only once and prior to calling other library functions; however, violations of this rule may be tolerated.

***cha\_version***

Returns a string that describes the current version of the CHAPRO library.

(char \*) **cha\_version**(void)

**Function arguments**

none

**Return Value**

Pointer to version string.

**Remarks**

An example of the return value, "CHAPro version 0.03, 6-Nov-2016".



***cha\_agc\_prepare***

Automatic-gain-control preparation function.

(int) **cha\_agc\_prepare**(CHA\_PTR **cp**, CHA\_DSL \***dsl**, CHA\_WDRC \***gha**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>dsl</b>	pointer to DSL prescription structure (see Appendix C)
<b>gha</b>	pointer to WDRC prescription structure (see Appendix D)

**Return Value**

Error code:  
0 – no error

**Remarks**

Initializes variables and allocates memory for automatic gain control. Chunk size is the number of samples read from the input signal and written to the output signal with each call to `cha_agc_process`.

**See Also**

**cha\_agc\_process**

***cha\_agc\_input***

Automatic-gain-control processing function.

(void) **cha\_agc\_input**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs single-channel, automatic-gain-control processing on CHAPRO input signal. Chunk size is the number of samples read from the input signal and written to the output signal.

**See Also**

**cha\_agc\_prepare**, **cha\_agc\_output**

***cha\_agc\_channel***

Automatic-gain-control processing function.

(void) **cha\_agc\_channel**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs multi-channel, automatic-gain-control processing. Chunk size is the number of samples read from the input signal and written to the output signal.

**See Also**

**cha\_agc\_prepare**

***cha\_agc\_output***

Automatic-gain-control processing function.

(void) **cha\_agc\_output**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs single-channel, automatic-gain-control processing on CHAPRO output signal. Chunk size is the number of samples read from the input signal and written to the output signal.

**See Also**

**cha\_agc\_prepare, cha\_agc\_input**

***cha\_compressor\_prepare***

Instantaneous-compression preparation function.

(int) **cha\_compressor\_prepare**(CHA\_PTR **cp**,  
float \***Lc**, float \***Gc**, double **lr**, int **np**, int **ds**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>Lc</b>	pointer to level array
<b>Lc</b>	pointer to gain array
<b>lr</b>	level reference
<b>np</b>	number of points in level and gain arrays
<b>ds</b>	down-sample factor

**Return Value**

Error code:  
0 – no error

**Remarks**

Initializes variables and allocates memory for instantaneous compression. Chunk size is the number of samples read from the input signal and written to the output signal with each call to `cha_compressor_process`.

**See Also**

**cha\_compressor\_process**

***cha\_compressor\_process***

Instantaneous-compression processing function.

(void) **cha\_compressor\_process**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs automatic-gain-control processing. Chunk size is the number of samples read from the input signal and written to the output signal.

**See Also**

**cha\_compressor\_prepare**

***cha\_afc\_prepare***

Save output signal for feedback management.

(int) **cha\_afc\_prepare**(CHA\_PTR **cp**,  
double **mu**, double **rho**, int **afl**, double **fbg**, int **sqm**)

**Function arguments**

<b>cp</b>	pointer to CHAPRO data structure
<b>mu</b>	AFC filter-estimation step size
<b>rho</b>	AFC filter-estimation forgetting factor
<b>fbg</b>	simulated-feedback gain between 0 and 1
<b>sqm</b>	option to save quality metric

**Return Value**

Error code:  
0 – no error

**Remarks**

Initializes variables and allocates memory for automatic gain control.

**See Also**

**cha\_afc\_input, cha\_afc\_output**

***cha\_afc\_input***

Process input signal to remove feedback.

(void) **cha\_afc\_input**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHAPRO data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Removes estimated feedback from input signal. Chunk size is the number of samples read from the input signal and written to the output signal. Optionally simulates feedback and saves misalignment error as a quality metric.

**See Also**

**cha\_afc\_prepare, cha\_afc\_output**



***cha\_afc\_output***

Save output signal for feedback management.

(void) **cha\_afc\_output**(CHA\_PTR **cp**, float \***x**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHAPRO data structure
<b>x</b>	pointer to input signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Assists feedback management by saving the output of the hearing-aid processing. Chunk size is the number of samples read from the input signal.

**See Also**

**cha\_afc\_prepare, cha\_afc\_input**

### ***cha\_cgtfb\_prepare***

Gammatone filter-bank preparation function.

(int) **cha\_cgtfb\_prepare**(CHA\_PTR **cp**, double \***fc**, double \***bw**,  
double **sr**, double **gd**, double **tw**, int **nc**, int **cs**)

#### **Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>fc</b>	pointer to list of center frequencies (Hz)
<b>bw</b>	pointer to list of bandwidths (Hz)
<b>sr</b>	sampling rate (samples/second)
<b>gd</b>	target group delay (ms)
<b>tw</b>	buffer length (ms) for determining zero gain
<b>nc</b>	number of channels
<b>cs</b>	chunk size

#### **Return Value**

Error code:  
0 – no error

#### **Remarks**

Initializes variables and allocates memory for the gammatone filter-bank. Chunk size is the number of samples read from the input and written to the output.

#### **See Also**

**cha\_cgtfb\_analyze**, **cha\_cgtfb\_synthesize**

***cha\_cgtfb\_analyze***

Gammatone filter-bank frequency-analysis function.

(void) **cha\_cgtfb\_analyze**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs automatic-gain-control processing. Chunk size is the number of samples read from the input and written to the output.

**See Also**

**cha\_cgtfb\_prepare**, **cha\_cgtfb\_synthesize**

***cha\_cgtfb\_synthesize***

Gammatone filter-bank frequency-synthesis function.

(void) **cha\_cgtfb\_synthesize**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs automatic-gain-control processing. Chunk size is the number of samples read from the input and written to the output.

**See Also**

**cha\_cgtfb\_prepare, cha\_cgtfb\_analyze**

***cha\_cfirfb\_prepare***

Complex FIR filter-bank preparation function.

(int) **cha\_cfirfb\_prepare**(CHA\_PTR **cp**, double \***cf**, int **nc**, double **sr**,  
int **nw**, int **wt**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>cf</b>	list frequency band edges (kHz)
<b>nc</b>	number of frequency bands
<b>sr</b>	sampling rate (samples/second)
<b>nw</b>	window size (samples)
<b>wt</b>	window type (0=Hamming, 1=Blackman)
<b>cs</b>	chunk size

**Return Value**

Error code:  
0 – no error

**Remarks**

Initializes variables and allocates memory for the complex FIR filter-bank. Chunk size is the number of samples read from the input and written to the output.

**See Also**

**cha\_cfirfb\_analyze**, **cha\_cfirfb\_synthesize**

***cha\_cfirfb\_analyze***

Complex FIR filter-bank frequency-analysis function.

(void) **cha\_cfirfb\_analyze**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to real input signal
<b>y</b>	pointer to complex output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs complex FIR filter-bank analysis. Chunk size is the number of samples read from the input and written to the output.

**See Also**

**cha\_cfirfb\_prepare**, **cha\_cfirfb\_synthesize**

***cha\_cfirfb\_synthesize***

Complex FIR filter-bank frequency-synthesis function.

(void) **cha\_cfirfb\_synthesize**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to complex input signal
<b>y</b>	pointer to real output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs complex FIR filter-bank synthesis. Chunk size is the number of samples read from the input and written to the output.

**See Also**

**cha\_cfirfb\_prepare, cha\_cfirfb\_analyze**

***cha\_firfb\_prepare***

FIR filter-bank preparation function.

(int) **cha\_firfb\_prepare**(CHA\_PTR **cp**,  
double \***cf**, int **nc**, double **sr**, int **nw**, int **wt**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>cf</b>	list frequency band edges (kHz)
<b>nc</b>	number of frequency bands
<b>sr</b>	sampling rate (samples/second)
<b>nw</b>	window size (samples)
<b>wt</b>	window type (0=Hamming, 1=Blackman)
<b>cs</b>	chunk size

**Return Value**

Error code:  
0 – no error

**Remarks**

Initializes variables and allocates memory for the FIR filter-bank. Chunk size is the number of samples read from the input and written to the output.

**See Also**

**cha\_firfb\_analyze, cha\_firfb\_synthesize**



***cha\_firfb\_analyze***

FIR filter-bank frequency-analysis function.

(void) **cha\_firfb\_analyze**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs FIR filter-bank analysis. Chunk size is the number of samples read from the input signal and written to the output signal.

**See Also**

**cha\_firfb\_prepare**, **cha\_firfb\_synthesize**

***cha\_firfb\_synthesize***

FIR filter-bank frequency-synthesis function.

(void) **cha\_firfb\_synthesize**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs FIR filter-bank synthesis. Chunk size is the number of samples read from the input signal and written to the output signal.

### ***cha\_iirfb\_prepare***

IIR filter-bank preparation function.

(int) **cha\_iirfb\_prepare**(CHA\_PTR **cp**,  
double \***z**, double \***p**, double \***g**, double \***d**, int **nc**, int **op**, double **sr**)

#### **Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>z</b>	pointer to IIR filter zeros
<b>b</b>	pointer to IIR filter poles
<b>g</b>	pointer IIR filter gain
<b>d</b>	pointer IIR filter delay
<b>nc</b>	number of frequency bands
<b>op</b>	number of filter coefficients for each band
<b>sr</b>	sampling rate (samples/second)

#### **Return Value**

Error code:  
0 – no error

#### **Remarks**

Initializes variables and allocates memory for the IIR filter-bank. The poles and zeros are complex numbers, so are stored as sequential real and imaginary parts. The size of the zeros & poles arrays (**z** & **p**) is two times the product of the number of frequency bands (**nc**) and the number of zeros per band (**nz**). The number of filter gains and filter delays (**g** & **d**) is the number of frequency bands (**nc**). Chunk size is the number of samples read from the input and written to the output.

#### **See Also**

**cha\_iirfb\_analyze**, **cha\_iirfb\_synthesize**

***cha\_iirfb\_analyze***

IIR filter-bank frequency-analysis function.

(void) **cha\_iirfb\_analyze**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs IIR filter-bank analysis. Chunk size is the number of samples read from the input signal and written to the output signal.

**See Also**

**cha\_iirfb\_prepare**, **cha\_iirfb\_synthesize**

***cha\_firfb\_synthesize***

IIR filter-bank frequency-synthesis function.

(void) **cha\_iirfb\_synthesize**(CHA\_PTR **cp**, float \***x**, float \***y**, int **cs**)

**Function arguments**

<b>cp</b>	pointer to CHA data structure
<b>x</b>	pointer to input signal
<b>y</b>	pointer to output signal
<b>cs</b>	chunk size

**Return Value**

none

**Remarks**

Performs IIR filter-bank synthesis. Chunk size is the number of samples read from the input signal and written to the output signal.

**See Also**

**cha\_iirfb\_prepare, cha\_iirfb\_analyze**

## Appendix A. Test programs

Several examples that test basic aspects of complex-gammatone filter-bank and instantaneous-compression.

- `tst_gfa` – tests complex-gammatone filter-bank analysis
- `tst_gfio` – tests simple waveform gammatone processing
- `tst_gfsc` – tests simple waveform gammatone processing with soundcard

Several examples that test basic aspects of FIR filter-bank and automatic-gain-control.

- `tst_ffa` – tests FIR filter-bank analysis
- `tst_ffio` – tests simple waveform FIR processing
- `tst_ffsc` – tests speech-waveform FIR & AGC with soundcard

Several examples that test basic aspects of IIR filter-bank with automatic-gain-control and adaptive feedback cancelation.

- `tst_ifa` – tests IIR filter-bank analysis
- `tst_ifio` – tests simple waveform IIR processing
- `tst_ifsc` – tests speech-waveform IIR & AGC
- `tst_iffb` – tests speech-waveform IIR & AFC with soundcard
- `tst_gha` – tests speech-waveform IIR & AFC & AGC with soundcard

These test programs are all written in C and produce results that are written to a subfolder called “test.” Each test program has a corresponding MATLAB script for viewing the results. The test programs that contain IIR processing require filter coefficients to be precomputed by a MATLAB script called `iirfb`.

## Appendix B. CLS Prescription

Structure CHA\_CLS specifies the CLS prescription.

```
#define CLS_MXCH 32          // maximum number of channels

typedef struct {
    int cm;                  // compression mode
    int nc;                  // number of channels
    double fc[CLS_MXCH];    // center frequency
    double bw[CLS_MXCH];    // bandwidth
    double Gcs[CLS_MXCH];   // gain at compression start
    double Gcm[CLS_MXCH];   // gain at compression middle
    double Gce[CLS_MXCH];   // gain at compression end
    double Gmx[CLS_MXCH];   // maximum gain
    double Lcs[CLS_MXCH];   // level at compression start
    double Lcm[CLS_MXCH];   // level at compression middle
    double Lce[CLS_MXCH];   // level at compression end
    double Lmx[CLS_MXCH];   // maximum output level
} CHA_CLS;
```

## Appendix C. DSL Prescription

Structure `CHA_DSL` specifies the DSL prescription.

```
#define DSL_MXCH 32                // maximum number of channels

typedef struct {
    double attack;                 // attack time (ms)
    double release;                // release time (ms)
    double maxdB;                  // maximum output (dB SPL)
    int ear;                       // 0=left, 1=right
    int nchannel;                  // number of channels
    double cross_freq[DSL_MXCH];  // cross frequencies (Hz)
    double tkgain[DSL_MXCH];      // compression-start gain
    double cr[DSL_MXCH];          // compression ratio
    double tk[DSL_MXCH];          // compression-start kneepoint
    double bolt[DSL_MXCH];        // broadband output limiting threshold
} CHA_DSL;
```



## Appendix D. WDRC Parameters

Structure `CHA_WDRC` specifies single-channel WDRC parameters

```
typedef struct {  
    double attack;           // attack time (ms)  
    double release;          // release time (ms)  
    double fs;               // sampling rate (Hz)  
    double maxdB;            // maximum signal (dB SPL)  
    double tkgain;           // compression-start gain  
    double tk;               // compression-start kneepoint  
    double cr;               // compression ratio  
    double bolt;             // broadband output limiting threshold  
} CHA_WDRC;
```